


8 KEYS DIFFERENTIAL COUNTER



- The 8 keys represent the normal schilling groups of leucocytes.
- Each key is the name and picture of particular cell.
- Bell rings at 100.
- The percentage of each type of cell is directly read from the window above each key.
- Total count available on the right window.
- Counterclockwise rotate handle will clear all figures.



Rue Fontaine St Pierre, 1C, Bat 6
B 5330 ASSESSE
BELGIUM

<http://www.amlab.be>

TEL : 32 83 66 89 30
FAX : 32 83 66 89 33

info@amlab.be